

Please amend the ABSTRACT OF THE DISCLOSURE to read as follows:

#### ABSTRACT OF THE DISCLOSURE

An electrosurgical biopsy device includes a stylet and a cannula movably mounted on a base. The stylet has a shaft with a head at its distal end and a stylet ablation element extending distally from the head. The stylet shaft is disposed through the cannula for axial translation therein between withdrawn and extended positions. The cannula has an opening at its distal end and a tissue cutting cannula ablation element adjacent the opening. Both ablation elements are activatable with energy that ablates adjacent tissue. A translation mechanism controllably moves (a) the stylet between the withdrawn and extended positions and (b) the cannula between a proximal position and a distal position relative to the base. ~~In use, with the stylet in the withdrawn position against the distal end of the cannula, and with the stylet ablation element activated, the stylet and the cannula are pushed through the skin and the underlying tissue until the stylet head is adjacent a targeted tissue mass. Next, the stylet is extended distally from the distal end of the cannula so that its head penetrates the tissue mass. The cannula ablation element is then activated, and the cannula is pushed through the tissue mass toward the stylet head, thereby cutting a "core" through the tissue mass that is captured as a tissue specimen within the distal end of the cannula. The cannula and the stylet are then removed from the patient's body.~~

**This Page is Inserted by IFW Indexing and Scanning  
Operations and is not part of the Official Record**

**BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ BLACK BORDERS
- ☐ IMAGE CUT OFF AT TOP, BOTTOM OR SIDES
- ☐ FADED TEXT OR DRAWING
- ☐ BLURRED OR ILLEGIBLE TEXT OR DRAWING
- ☐ SKEWED/SLANTED IMAGES
- ☐ COLOR OR BLACK AND WHITE PHOTOGRAPHS
- ☐ GRAY SCALE DOCUMENTS
- ☒ LINES OR MARKS ON ORIGINAL DOCUMENT
- ☐ REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY
- ☐ OTHER: \_\_\_\_\_

**IMAGES ARE BEST AVAILABLE COPY.**

**As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.**